

LUXVISION TABLES CATALOG

ELECTRICAL TABLES FOR ONE INSTRUMENT



FUNCTIONALITY AND ENDURANCE IN ONE SINGLE TABLE





LUXVISION ET-SERIES

ELECTRICAL TABLES FOR ONE INSTRUMENT

Functionality and endurance
in one single table.

ET-150

Electrical Table



* Equipment not included

ET-150

Electrical Table

The Luxvision ET-150 is ideal for eye care practice and easily accommodates one instrument. The streamlined design offers a movable base with locking wheels and offers height adjustment. High-grade casters ensure safe and stable movement. The plug for the switch of raising and landing under the table.

Features

- Convenient Equipment Accommodation
- Locking Wheels
- Adjustable Height
- Quality Materials and Design

Convenient Equipment Accommodation

The ET-150 Table allows you to conveniently operate almost any ophthalmic instrument. The table features an AC connection and power switch.

ET-160

Electrical Table

ET-160

Electrical Table

The Luxvision ET-160 is ideal for eye care practice and easily accommodates one instrument. The streamlined design offers a moveable base with locking wheels and offers height adjustment. High-grade casters ensure safe and stable movement. The plug for the switch of raising and landing under the table.

Features


- Convenient Equipment Accommodation
- Locking Wheels
- Adjustable Height
- Quality Materials and Design



* Equipment not included

ET-175

Electrical Table



* Equipment not included

ET-175

Electrical Table

The Luxvision ET-175 is ideal for eye care practice and easily accommodates one instrument. The streamlined design offers a movable base with locking wheels and offers motorized height adjustment controls on the pedestal. High-grade casters ensure safe and stable movement.

Features

- Motorized Height
- Locking Wheels
- Convenient Equipment Accommodation
- Quality Materials and Design

OPTIONS FOR 1 and 2 EQUIPMENTS COMPATIBLE WITH THE ET-175

*Top table only. Does not include equipment



GET-15

Short Table Top

Features

The Short TBL is adaptable to the ET-175 and ET-185 tables. It can hold one equipment.

- 20.87 in x 14.17 in (530 mm x 360 mm)
- Quality Materials and Design



GET-15

Long Table Top

Features

The Long TBL is adaptable to the ET-175 and ET-185 tables. It can hold up to two equipments.

- 32 in x 18 in (812.8 mm x 457.2 mm)
- Quality Materials and Design

ET-700

Electrical Table



*Equipment not included

ET-700

Electrical Table

The Luxvision ET-700 is ideal for eye care practice and easily accommodates one instrument. It is also designed with a sliding base to improve accessibility. Locking wheels and stable height adjustment make this electric table an excellent choice. It can be used for disabled patients.

SUITABLE FOR PATIENTS WITH DISABILITIES

*Table does not include equipment



Features

• Adjustable Height

Patients won't have to reach, bend or stretch to meet the instrument. The ET-700 offers an excellent range of height adjustment. It is also wheelchair accessible.

• Convenient Equipment Accommodation

The single instrument table ET-700 allows you to conveniently operate almost any ophthalmic instrument. The table features an AC connection and power switch.

• Locking Wheels

The table provides easy access and is wheelchair accessible. Simply roll the table to your desired location and lock the wheels as needed.

• Quality Materials and Design

The ET-700 was crafted using quality materials, including high-grade casters to ensure safe and stable movement.

Motorized Height Adjustment

Patients won't have to bend or stretch to meet the instrument. The ET-700 offers an excellent range of motorized height adjustment. Controls are located on the pedestal.

SPECIFICATIONS

ET-150



Table Size	23.62 - 15.75 in (600 - 400 mm)
Table Height of Raise and Landing	24.80 in - 32.28 in (630 - 820 mm)
Input Voltage	220V/50Hz - 110V/60Hz
Max Loading	121 lb (55 kg)
Weight	48.40 lb (22 kg)

ET-160



Table Size	15.75 x 22.83 in (400 X 580 mm)
Table Height of Raise and Landing	24.80 in - 32.28 in (630 - 820 mm)
Input Voltage	220V/50Hz - 110V/60Hz
Max Loading	165 lb (75 kg)
Weight	42.90 lb (19.5 kg)

ET-175



Dimensions	20.07 (W) x 17.72 (D) in (510 (W) x 450 (D) mm)
Table minimum height/stroke	25.99 in / 6.30 in (660 mm/160 mm)
Weight	44 lb (20 Kg (body only))
Maximum load	88 lb (40 Kg)
Input power supply (220V system)	AC220V, 230V, 240V
Input power supply (110V system)	AC110V, 120V
Input power frequency	50 / 60Hz
Power consumption	220VA
Electric safe standard	Conform to Standard 1 EC601-1 Class 1 Type B

ET-700



Table Size	24.80 in x 18.11 in (630 mm x 460 mm)
Range of Table Elevation	24.80 in x 32.28 in (630 mm - 820 mm)
Max Loading	121 lb (55 Kg)
Power	220V or 110V/60Hz