

LCP-7800

CHART PROJECTOR

TO PERFORM PRECISE ACUITY TEST WITH HIGH QUALITY CHARTS, LUXVISION CHART PROJECTOR CAN BE THE RIGHT ANSWER. LCP-7800 delivers its promises by utilizing the total refraction system to perform any test you want accurately and efficiently. The various kinds of chart can be changed more precisely within very short time.

SPECIFICATIONS

Chart	41 chart 34 mask Red/Green & Polarization Filters
Projection Distance	2.5-8m
Chart Rotation Speed Average	0.15 sec
Projection Magnification	30 x at 5m
Power Saving Program	Automatic Switch off (10min) 2 programs with a maximum of 30 charts each
Tilt Angle	15°
Power Supply	110-220V / 200-240V, 50 / 60 Hz
Power Consumption	Max 0.6A
Lamp	LED 4W
Dimension	270(W) x 182(D) x 230(H) mm
Weight	3.5kg



LCP-7800

CHART PROJECTOR

LCP-7800

REMARKABLE FEATURES



CLEAR AND FAST HIGH QUALITY CHARTS

LCP-7800 has high resolution power of line, so provides neat and clean graphics and a fast chart changing speed making it easy for operators to use and comfortable for customers.

EASY AND CONVENIENT/ BEST NETWORKING

Any refraction system including Luxvision digital refractor can incorporate LCP-7800 to provide fast and efficient vision test.

BRIGHT AND EFFICIENT LED LIGHT

Semi-permanent lifetime LED light source can provide brighter and clearer charts than conventional bulb light.

VARIOUS AND PERFECT WIDE RANGE OF FEATURES

LCP-7800 with 41 charts including Red/Green and polarizing filters are the new solution for Binocular balance test, Stereo test, Anisikonia test and Fusion test.